

Datasheet: NEB-Cterm

Description: Rabbit polyclonal Antibody (batch#6964)

Specificity: Nebulin C-terminal 35kDa region

Other names: NEB-C rab ab
Product Type: Polyclonal Antibody
Isotype: Polyclonal IgG
Quantity: 100µL (70 µg)

Product Details

This antiserum has been raised to the carboxy-terminal 35kDa region of nebulin. This part of nebulin is not composed of highly repetitive sequences (1,2). It therefore provides a suitable region for specific epitope recognition. Structurally, this epitope localizes inside Z-disks within the barbed-end region of the thin filament (3). Functionally, this segment of nebulin has been implicated in binding to titin (4), to myopalladin (5), and to N-WASP after stimulation with IGF-1 (see (6) and related product: Nebulin Cterm – phosphospecific antibody).

The NEB-Cterm antibody has been affinity-purified with the specific antigen coupled to an affinity resign. Affinity-purified IgGs after elution are provided at a concentration of 0.7 μ g/ μ l. For long-term storage, aliquoting, snap-freezing and storage at -80°C is recommended. For storage for up to 6 months, storage at 4°C without re-freezing is recommended.

The specificity of this antibody has been verified by Western blots by comparing wildtype and Nebulin-KO mice: Inactivation of the nebulin gene in mice by gene targeting leads to nebulin-deficient myofibrils (7 8). Western blots with NEB-Cterm antibody confirm lack of nebulin expression in NEB-KO mice (Figure 1, right).


In human, nebulin-deficiency is causative for the mainstream form nebulin myopathy (9). Therefore, nebulin Cterm antibodies are useful to detect nebulin deficiency (or truncations (10-12) in patients with nemaline myopathies.

Rabbit IgG polyclonal NEB-Cterm:

Recommended usages: Western blots, immunofluorescence

Tested species: Human and mouse

Related Products: Nebulin-Nterm antibody, Nebulin-Cterm antibody (phospho-specific)



WT KO WT KO WT KO

Figure 1. Western blot characaterization of anti-NEB-Cterm antibodies.

Simultaneous detection of titin and nebulin in the two channel-mode using different titin and nebulin antibodies (see also 8). Anti-NEB-Cterm detects nebulin in wild-type but not nebulin KO mice (Figure kindly provided by Danielle Buck and Henk Granzier, University of Arizona)




Figure 2.

Immunofluorescence labeling of isolated myofibrils with NEB-Cterm antibody (top; green).

Middle: Labeling of sarcomeric Z-disks with desmin (center, in red). The merge of both channels assigns the staining of NEB-Cterm to the sarcomeric Z-disk (bottom: merge in yellow; courtesy of Akira Hanashima, Chiba University, Japan).

Literature

- Labeit S, Kolmerer B. The complete primary structure of human nebulin and its correlation to muscle structure. J Mol Biol. 1995;248:308-15.
- 2. Labeit S, Ottenheijm CA, Granzier H. Nebulin, a major player in muscle health and disease. FASEB J. 2011;25:822-9.
- 3. Millevoi S, Trombitas K, Kolmerer B, Kostin S, Schaper J, Pelin K, Granzier H, Labeit S. Characterization of nebulette and nebulin and emerging concepts of their roles for vertebrate Z-discs. J Mol Biol. 1998;282:111-123.
- Labeit S, Lahmers S, Burkart C, Fong C, McNabb M, Witt S, Witt C, Labeit D, Granzier H. Expression of distinct classes of titin isoforms in striated and smooth muscles by alternative splicing, and their conserved interaction with filamins. J Mol Biol. 2006; 362: 664-681.
- Bang ML, Mudry RE, McElhinny AS, Trombitás K, Geach AJ, Yamasaki R, Sorimachi H, Granzier H, Gregorio CC, Labeit S. Myopalladin, a novel 145-kilodalton sarcomeric protein with multiple roles in Z-disc and I-band protein assemblies. J Cell Biol. 2001; 153: 413-427.
- 6. Takano K, Watanabe-Takano H, Suetsugu S, Kurita S, Tsujita K, Kimura S, Karatsu T, Takenawa T, Endo T. Nebulin and N-WASP cooperate to cause IGF-1-induced sarcomeric actin filament formation. Science. 2010; 330:1536-1540.
- 7. Bang ML, Li X, Littlefield R, Bremner S, Thor A, Knowlton KU, Lieber RL, Chen J. Nebulin-deficient mice exhibit shorter thin filament lengths and reduced contractile function in skeletal muscle. J Cell Biol. 2006;173:905-16.
- 8. Witt CC, Burkart C, Labeit D, McNabb M, Wu Y, Granzier H, Labeit S. Nebulin regulates thin filament length, contractility, and Z-disk structure in vivo. EMBO J. 2006;25:3843-55.
- Pelin K, Hilpelä P, Donner K, Sewry C, Akkari PA, Wilton SD, Wattanasirichaigoon D, Bang ML, Centner T, Hanefeld F, Odent S, Fardeau M, Urtizberea JA, Muntoni F, Dubowitz V, Beggs AH, Laing NG, Labeit S, de la Chapelle A, Wallgren-Pettersson C. Mutations in the nebulin gene associated with autosomal recessive nemaline myopathy. Proc Natl Acad Sci U S A. 1999;96:2305-10.
- 10. Ottenheijm CA, Witt CC, Stienen GJ, Labeit S, Beggs AH, Granzier H. Thin filament length dysregulation contributes to muscle weakness in nemaline myopathy patients with nebulin deficiency. Hum Mol Genet. 2009;18:2359-69.
- 11. Lawlor MW, Ottenheijm ČA, Lehtokari VL, Cho K, Pelin K, Wallgren-Pettersson C, Granzier H, Beggs AH. Novel mutations in NEB cause abnormal nebulin expression and markedly impaired muscle force generation in severe nemaline myopathy. Skelet Muscle. 2011;1:23.
- 12. Sewry CA, Brown SC, Pelin K, Jungbluth H, Wallgren-Pettersson C, Labeit S, Manzur A, Muntoni F. Abnormalities in the expression of nebulin in chromosome-2 linked nemaline myopathy. Neuromuscul Disord. 2001;11:146-53.

For research purposes only, unless otherwise specified in writing by Myomedix GmbH.



Myomedix GmbH • MIttelstr. 42 • 68169 Mannheim

Tel.: +49621 3832716 • Email: info@myomedix.com

Web: www.myomedix.com

Printed on 29.10.2012